Nanopore sequencing

Nick Loman

Introduction to nanopore

- MinION is the size of a:
 - USB stick
 - "around the same length as an iPhone" [1]
 - "no larger than a typical smartphone" [2]
 - Fits in your pocket [3]
 - "10 cm \times 2 cm \times 3.3 cm; approximately 90 g" [4]
 - Office stapler [5]
 - Ruler
 - A 4" ruler [6]
 - Chocolate bar [7]



Size comparison of one of the @nanopore MinIONs and a banana. Actually, search youtube for 'minion and banana'. :-)



[1] Loman et al. Bioinformatics 2014, [2] Quick et al. GigaScience 2014, [3] Urban et al. bioRxiv 2015, [4] Wang et al. Front. Microbiol 2015 [5] Risse, Watson et al. bioRxiv 2015, [6] Unpublished, probably forever, [7] Manuscript under preparation

Priming

- Nanopore sequencing starts when:
 - Fuel is added to flow cell
 - ATP makes ratcheting enzyme work
 - More ATP faster speed, more data
 - Currently 35bp/s, due to change to 70bp/s soon
 - Potentially 500 bp/s or higher?
 - Electrical potential applied
 - Movement of ions
 - "lonic battery"
 - Controlled by Python scripts
 - Bias voltage changes over time but quality drops off

QC

- MinION Mk I flow cell is actually
 - 2000 wells, theoretically each containing 1 pore
 - But Poisson loading means some wells >1 pore (unusable – "multiple") or 0 pores ("None")
 - 1 sensor for each 4 well cluster can address one
 - Each sensor is a miniature axopatch
 - Pores divided into groups
 - "Remuxing" after length of time
 - Again controlled by Python

Interrogating MinKNOW

- Read length histogram
 - Counts 'events' events ~ 1base (but in practice <1 due to skips)</p>
 - If read is 2D then 2x events- don't get too excited
- Pore view

Poretools

Progress monitoring

Metrichor

Real-time analysis

What's in my pot?

Extract FASTA files

Hybrid assembly with SPAdes

Compare with Illumina only assembly

Nanopore-only assembly

- Is possible!
 - Nanocorrect
 - Celera
 - Nanopolish
 - Resulting assembly is 99.5% accurate at 30x mainly homopolymer errors
 - Rather time consuming at present

QUAST again

Inspecting alignments with BWA and IGV

Do it again